



DIGITAL IMAGING MADE EASY

## **QCapture Pro 7™ Information**

QCapture Pro 7 software is available for download and will run as a 30 day free trial. After 30 days the software will revert to demonstration mode.

QImaging recommends an evaluation of the software be completed within two weeks and if results are satisfactory, the order to purchase be submitted. This ensures the security dongle will arrive before the trial expires, enabling continued use of the software in full mode without disruption.

### **Details:**

Part Number: 01-QCAP-PRO7-S

Part Description: QCAPTURE PRO7 CD W/DONGLE

Price: \$250.00 USD each

### **Addition Information:**

- Download [QCapture Pro 7](#) software and complete the installation.
- Installation of the [latest QCam camera driver](#) is required.
- View the [Basic Concepts of QCapture Pro 7 Software](#) video on QITube.
- When installing on a new computer, the PCI FireWire card must be transferred from the old PC to the new PC. Specific third party PCI FireWire cards are tested and supported. The PCI FireWire cards that QImaging currently ships with new cameras include:
  - [PCIe FireWire Card](#)
  - [PCI FireWire Card](#)
- QCapture Pro 7 is not supported on Windows 7 – Enterprise. For additional information, refer to the [QCapture Pro 7 Datasheet](#).

### **To Order:**

To order QCapture Pro 7 software, [complete and submit the online request form](#).

19535 56<sup>th</sup> Avenue, Suite 101

Surrey, BC

Canada V3S 6K3

tel 604.530.5800

fax 604.648.8277

[www.qimaging.com](http://www.qimaging.com)